



| Description | Code | Masse/Weight Kg/p | Ep./Th. mm | l./w. mm |
|--------------------|--------|----------------------|---------------|-------------|
| CO SLT NERV 50 I3 | 425189 | 0,163 | 0,80 | 57 |
| CO SLT NERV 75 I3 | 425194 | 0,174 | 0,80 | 81 |
| CO SLT NERV 100 I3 | 425199 | 0,196 | 0,80 | 105 |
| CO SLT NERV 123 I3 | 425204 | 0,250 | 0,80 | 129 |
| CO SLT NERV 147 I3 | 425209 | 0,294 | 0,80 | 153 |
| CO SLT NERV 195 I3 | 425219 | 0,304 | 0,80 | 201 |
| CO SLT NERV 316 I3 | 425234 | 0,424 | 0,80 | 321 |
| CO SLT NERV 412 I3 | 425244 | 0,511 | 0,80 | 417 |
| CO SLT NERV 507 I3 | 425264 | 0,620 | 0,80 | 513 |
| CO SLT NERV 604 I3 | 425274 | 0,717 | 0,80 | 609 |

SS316L : Stainless steel according to NF EN 10088-2 (316L).
I3 : Tôle d'acier inoxydable suivant NF EN 10088-2 (316L).



www.niedaxfrance.fr

COUVERCLE DERIVATION I3
 TAKE OFF COVER SS316L



ASSISTANCE TECHNIQUE
 +33 (0)1 55 95 92 17
 pole.technique@niedax.fr

FT02492

NIEDAX
 France

